

Sulfachloropyridazine

Catalog No: #S2184


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Sulfachloropyridazine
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	284.72
Formulation	C10H9ClN4O2S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
 SMILES: NC1=CC=C(C=C1)S(=O)(=O)NC1=NN=C(Cl)C=C1
 Pathways: Others
 Receptor: Others
 Boiling pt: 559.7°C
 Melting pt: 186-187°C
 Solubility: Soluble in DMSO
 Appearance: Yellow Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Alavi FK, et al. J Vet Pharmacol Ther. 1993 Jun;16(2):232-6.

Note: This product is for in vitro research use only