

Sulfamethazine sodium salt

Catalog No: #S2184L


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Sulfamethazine sodium salt
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	300.31
Formulation	C ₁₂ H ₁₃ N ₄ NaO ₂ S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology SMILES: [Na+].CC1=CC(C)=NC([N-]S(=O)(=O)C2=CC=C(N)C=C2)=N1 Pathways: Others Receptor: Others Boiling pt: 526.2°C Melting pt: 176°C Solubility: Soluble in Water Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Fukahori S, et al. J Environ Manage. 2015 Jul 1;157:103-10.

Note: This product is for in vitro research use only