

Sodium L-ascorbyl-2-phosphate

Catalog No: #S2185


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Sodium L-ascorbyl-2-phosphate
Brief Description	Inhibitors
Purification	98.00%
Target Name	Vitamin inhibitor
Calculated MW	322.05
Formulation	C6H6Na3O9P
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Metabolism

systemSMILES:C([C@@H])([C@@H]1C(=C(C(=O)O1)OP(=O)([O-])[O-])O)O.[Na+].[Na+]Pathways:MetabolismReceptor:VCBoiling pt:Melting pt:84CSolubility:Soluble in WaterAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Ares AM, et al. J Chromatogr A. 2013 Oct 25;1313:78-95.

Note: This product is for in vitro research use only