

Tris(2,2,2-trifluoroethyl) borate

Catalog No: #S2186


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Tris(2,2,2-trifluoroethyl) borate
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	307.91
Formulation	C6H6BF9O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: FC(F)(F)COB(OCC(F)(F)F)OCC(F)(F)F
Pathways: Others
Receptor: Others
Boiling pt: 94.3°C
Melting pt: Solubility: DMSO: > 11 mM
Appearance: Liquid
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Rait JL. Aust N Z J Ophthalmol. 1999, 27(1):57-64.

Note: This product is for in vitro research use only