

Tacrine hydrochloride hydrate

Catalog No: #S2191


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Tacrine hydrochloride hydrate
Brief Description	Inhibitors
Purification	99.89%
Target Name	AChR inhibitor
Calculated MW	234.72
Formulation	C ₁₃ H ₁₄ N ₂ HCl·xH ₂ O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
 SMILES: O=C1NC1=C2CCCCC2=NC2=CC=CC=C2
 Pathways: Neuroscience
 Receptor: AChR
 Boiling pt: Melting pt: Solubility: Soluble in Water
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Wang H, et al. Zhongguo Yao Li Xue Bao. 1998, 19(1):27-30.

Note: This product is for in vitro research use only