

Indazole


Catalog No: #S2192

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Indazole
Brief Description	Inhibitors
Purification	98.00%
Target Name	Sodium Channel inhibitor
Calculated MW	118.14
Formulation	C7H6N2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: N1N=CC2=CC=CC=C12
Pathways: Membrane transporter/Ion channel
Receptor: Sodium Channel
Boiling pt: 270°C
Melting pt: 146-150°C
Solubility: DMSO: > 10 mM
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.
Stock solution can be stored below -20 °C for several month.

References

1. Matsumura N, et al. *Epilepsy Res.* 2013 May;104(3):203-16.

Note: This product is for in vitro research use only