

2-Chloropyrazine


Catalog No: #S2193

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	2-Chloropyrazine
Brief Description	Inhibitors
Purification	98.00%
Target Name	CDK inhibitor
Calculated MW	114.53
Formulation	C4H3ClN2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: C1=NC=CN=C1Cl
Pathways: Cell Cycle/Checkpoint Receptor: CDK2
Boiling pt: 155.3°C
Melting pt: Solubility: DMSO: 10 mM
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.
Stock solution can be stored below -20 °C for several months.

References

1. Bgouin JM, Gosmini C. J Org Chem. 2009 Apr 17;74(8):3221-4.

Note: This product is for in vitro research use only