

Epiandrosterone

Catalog No: #S2195


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-------------------|--|
| Product Name | Epiandrosterone |
| Brief Description | Inhibitors |
| Purification | 99.00% |
| Target Name | Androgen Receptor agonist; Estrogen/progestogen Receptor |
| Calculated MW | 290.44 |
| Formulation | C19H30O2 |
| Storage | 3 years -20°C powder; 2 years -80°C in solvent; |

Images



Product Description

Research Area:Endocrine

systemSMILES:C[C@]12CC[C@@H](C[C@@H]1CC[C@@H]3[C@@H]2CC[C@]4([C@H]3CCC4=O)C)OPathways:Endocrinology/HormonesReceptor:Androgen Receptor; PRBoiling pt:413.1CMelting pt:172-174CSolubility:Ethanol: 58 mg/mL (199.69 mM); DMSO: 29 mg/mL warmed (99.84 mM)Appearance:Gray or White Crystal PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Gupte SA, et al. J Pharmacol Exp Ther, 2002, 301(1), 299-305.

Note: This product is for in vitro research use only