

## 10-Deacetylbaocatin III

Catalog No: #S2197


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	10-Deacetylbaocatin III
Brief Description	Inhibitors
Purification	99.50%
Target Name	Others
Calculated MW	544.59
Formulation	C29H36O10
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:OtherSMILES:CC1=C2[C@H](C(=O)[C@@]3([C@H](C[C@@H]4[C@]([C@H]3[C@@H]([C@@](C2(C)C)(C[C@@H]1O)O)OC(=O)C5=CC=CC=C5)(CO4)OC(=O)C)O)C)OPathways:OthersReceptor:OthersBoiling pt:716.8CMelting pt:231-236CSolubility:DMSO: 109 mg/mL (200.15 mM)Appearance:White SolidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Zhang JF, et al. Cell Biol Int Rep (2010). 2010; 18(1):e00009.

Note: This product is for in vitro research use only