

## Myclobutanil


Catalog No: #S2199

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Myclobutanil
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others inhibitor
Calculated MW	288.78
Formulation	C15H17ClN4
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:ImmunologySMILES:CCCC(CN1C=NC=N1)(C#N)C1=CC=C(Cl)C=C1Pathways:OthersReceptor:OthersBoiling pt:202-208CMelting pt:63-68CSolubility:Soluble in WaterAppearance:Faint Yellow SolidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Yan J, et al. Chirality. 2014 Jan;26(1):51-5.

Note: This product is for in vitro research use only