

L-Carnitine hydrochloride


Catalog No: #S2203

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	L-Carnitine hydrochloride
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	197.66
Formulation	C7H15NO3HCl
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Metabolism systemSMILES:C[N+](C)(C)C[C@H](CC(=O)O)O.[Cl-]

Pathways:OthersReceptor:CRATBoiling pt:Melting pt:Solubility:DMSO: > 10 mMAppearance:White CrystalPowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Merra G, et al. World J Gastroenterol. 2012 Sep 28;18(36):5065-71.

Note: This product is for in vitro research use only