

## Diaveridine

Catalog No: #S2204


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Diaveridine
Brief Description	Inhibitors
Purification	99.19%
Target Name	Others
Calculated MW	260.29
Formulation	C13H16N4O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Immunology  
SMILES: COC1=C(OC)C=C(CC2=CN=C(N)N=C2N)C=C1  
Pathways: Others  
Receptor: Others  
Boiling pt: 506.1°C  
Melting pt: 232 - 236°C  
Solubility: DMSO: > 10 mM  
Appearance: White Solid  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Ono T, et al. Environ Toxicol Pharmacol. 1997 Sep;3(4):297-306.

Note: This product is for in vitro research use only