

Edetol

Catalog No: #S2206


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Edetol
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	292.41
Formulation	C ₁₄ H ₃₂ N ₂ O ₄
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: CC(O)CN(CC(C)O)CC(C)O
Pathways: Others
Receptor: Others
Boiling pt: 175-181°C
Melting pt: 32°C
Solubility: Soluble in Water
Appearance: Transparent Colorless - Faint Yellow Liquid
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

N/A

Note: This product is for in vitro research use only