

Adrenosterone

Catalog No: #S2207


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Adrenosterone
Brief Description	Inhibitors
Purification	99.00%
Target Name	Androgen Receptor antagonist
Calculated MW	300.4
Formulation	C19H24O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Endocrine

systemSMILES:C[C@]12CCC(=O)C=C1CC[C@@H]3[C@@H]2C(=O)C[C@]4([C@H]3CCC4=O)CPathways:Endocrinology/HormonesReceptor:Androgen ReceptorBoiling pt:476.1CMelting pt:219-222CSolubility:DMSO: 10 mMAppearance:Gray or Faint Yellow SolidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Blasco M, et al. J Chromatogr B Analyt Technol Biomed Life Sci. 2009 May 15;877(14-15):1509-15.

Note: This product is for in vitro research use only