

## Geniposide

Catalog No: #S2212


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Geniposide
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	388.37
Formulation	C17H24O10
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Nervous

systemSMILES:COC(=O)C1=CO[C@H]([C@H]2[C@@H]1CC=C2CO)O[C@H]3[C@@H]([C@H]([C@@H]([C@H](O3)CO)O)O)OPathways:OthersReceptor:OthersBoiling pt:641.4CMelting pt:161-162CSolubility:Water: 78 mg/mL (200.83 mM); DMSO: 78 mg/mL (200.83 mM)Appearance:White CrystalRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Zhang Y, et al. Neurochem Int. 2015 Apr 13. pii: S0197-0186(15)00063-7.

Note: This product is for in vitro research use only