

Terconazole

Catalog No: #S2213


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Terconazole
Brief Description	Inhibitors
Purification	99.94%
Target Name	Antifection agonist
Calculated MW	532.46
Formulation	C26H31Cl2N5O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: CC(C)N1CCN(CC1)C2=CC=C(C=C2)OC[C@H]3CO[C@](O3)(CN4C=NC=N4)C5=C(C=C(C=C5)Cl)Cl Pathways: Microbiology & virology Receptor: 14- α Demethylase Boiling pt: Melting pt: Solubility: Soluble in DMSO Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Zhao XY, et al. Chin Med J (Engl). 2012 Mar;125(6):1068-75.

Note: This product is for in vitro research use only