

Flumethrin

Catalog No: #S2214


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Flumethrin
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	510.38
Formulation	C28H22Cl2FNO3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology
 SMILES: CC1(C(C1C(=O)OC(C#N)C2=CC(=C(C=C2)F)OC3=CC=CC=C3)/C=C(/C4=CC=C(C=C4)Cl)\Cl)Cl
 Pathways: Others
 Receptor: Antiparasitic
 Boiling pt: Melting pt: Solubility: Soluble in Water
 Appearance: Liquid
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Dantas-Torres F, et al. Parasit Vectors. 2013 Aug 23;6(1):245.

Note: This product is for in vitro research use only