

Cephalomannine

Catalog No: #S2217


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cephalomannine
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	831.91
Formulation	C45H53NO14
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:CancerSMILES:C/C=C(\C)/C(=O)N[C@@H](C1=CC=CC=C1)[C@H](C(=O)O[C@H]2C[C@]3([C@H]([C@H]4[C@@]([C@H](C[C@@H]5[C@]4(CO5)OC(=O)C)O)(C(=O)[C@@H](C=C2C)C3(C)C)OC(=O)C)OC(=O)C6=CC=CC=C6)O)OPathways:OthersReceptor:OthersBoiling pt:929.5CMelting pt:139-142CSolubility:Ethanol: 100 mg/mL (120.2 mM); DMSO: 100 mg/mL (120.2 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Helson L. Int J Oncol. 1993 Feb;2(2):297-9.

Note: This product is for in vitro research use only