

Kinetin


Catalog No: #S2223

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Kinetin
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	215.21
Formulation	C10H9N5O
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:C(NC1=NC=NC2=C1NC=N2)C1=CC=CO1Pathways:OthersReceptor:OthersBoiling pt:524.4CMelting pt:269-271CSolubility:DMSO: 8 mg/mL warmed (37.17 mM)Appearance:Gray Amorphous PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Chiu PC, et al. J Cosmet Dermatol, 2007, 6(4), 243-249.

Note: This product is for in vitro research use only