

4-Hydroxyantipyrine

Catalog No: #S2225


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	4-Hydroxyantipyrine
Brief Description	Inhibitors
Purification	99.00%
Target Name	COX
Calculated MW	204.23
Formulation	C11H12N2O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Nervous systemSMILES:CN1N(C(=O)C(O)=C1C)C1=CC=CC=C1Pathways:Immunology/Inflammation; NeuroscienceReceptor:COX-1; COX-2Boiling pt:314.3CMelting pt:184-186CSolubility:DMSO: 10 mMAppearance:Yellow PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Brune K, Neubert A. Clin Exp Rheumatol. 2001 Nov-Dec;19(6 Suppl 25):S51-7.

Note: This product is for in vitro research use only