

Pergolide mesylate salt

Catalog No: #S2226


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Pergolide mesylate salt
Brief Description	Inhibitors
Purification	98.00%
Target Name	5-HT Receptor antagonist; Adrenergic Receptor agonist; Dopamine Receptor agonist
Calculated MW	410.59
Formulation	C20H30N2O3S2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Others SMILES: CCCN1C[C@@H](C[C@H]2[C@H]1CC3=CNC4=CC=CC2=C34)CSC.CS(=O)(=O)O Pathways: Neuroscience Receptor: 5-HT

receptor; Adrenergic receptor; Dopamine receptor Boiling pt: Melting pt: 252-254°C Solubility: DMSO 10 mg/mL; Water 82 mg/mL Appearance: Off-white

Color Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Fuller RW, et al. J Pharm Pharmacol. 1985 Apr;37(4):268-70.

Note: This product is for in vitro research use only