

## Hemicholinium-3

Catalog No: #S2227


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Hemicholinium-3
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others inhibitor
Calculated MW	574.35
Formulation	C <sub>24</sub> H <sub>34</sub> BrN <sub>2</sub> O <sub>4</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area:Nervous

systemSMILES:[Br-].[Br-].C[N+](C)CCOC(O)(C1)C1=CC=C(C=C1)C1=CC=C(C=C1)C1(O)C[N+](C)(C)CCO1Pathways:OthersReceptor:ChRBoling  
 pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Soukup O, et al. Toxicol Lett. 2012 Aug 3;212(3):315-9.

Note: This product is for in vitro research use only