

## Adipic acid

Catalog No: #S2228


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Adipic acid
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	146.14
Formulation	C6H10O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:OtherSMILES:OC(=O)CCCCC(=O)OPathways:OthersReceptor:OthersBoiling pt:338.5CMelting pt:151-153CSolubility:DMSO: 10 mMAppearance:White CrystalRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Han L, et al. Sheng Wu Gong Cheng Xue Bao. 2013 Oct;29(10):1374-85.

Note: This product is for in vitro research use only