

17-Beta-Estradiol-3,17-Dipropionate

Catalog No: #S2231


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	17-Beta-Estradiol-3,17-Dipropionate
Brief Description	Inhibitors
Purification	98.89%
Target Name	Estrogen/progestogen Receptor agonist
Calculated MW	384.52
Formulation	C24H32O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:CancerSMILES:CCC(=O)O[C@H]1CC[C@@H]2[C@@]1(CC[C@H]3[C@H]2CCC4=C3C=CC(=C4)OC(=O)CC)CPathways:Endocrinology/HormonesReceptor:ERBoling pt:462.4CMelting pt:104-106CSolubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Noguchi M, et al.
J Vet Med Sci. 2013; 75(3):343-8.

Note: This product is for in vitro research use only