

## Cariporide

Catalog No: #S2238


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Cariporide
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others inhibitor
Calculated MW	283.35
Formulation	C12H17N3O3S
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Cardiovascular SystemSMILES:CC(C)C1=C(C=C(C=C1)C(=O)N=C(N)N)S(C)(=O)=OPathways:OthersReceptor:NHE2Boiling pt:542.8CMelting pt:Solubility:Soluble in DMSOAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Yang X, et al. J Surg Res. 2013 Dec;185(2):797-804.
2. Wu S, et al. Pharmacology. 2013; 91(3-4):123-30.

Note: This product is for in vitro research use only