


Description

| | |
|-------------------|--|
| Product Name | DY-131 |
| Brief Description | Inhibitors |
| Purification | 98.00% |
| Target Name | Estrogen/progestogen Receptor agonist |
| Calculated MW | 311.38 |
| Formulation | C18H21N3O2 |
| Storage | 3 years -20°C powder;2 years -80°C in solvent; |

Images



Product Description

Research Area:Endocrine

systemSMILES:CCN(CC)C1=CC=C(C=C1)/C=N/NC(=O)C2=CC=C(C=C2)OPathways:Endocrinology/HormonesReceptor:ERRβ; ERRγBoling

pt:Melting pt:Solubility:DMSO: 14 mg/mLAppearance:Yellow PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1.Wang Y, et al. ACS Chem Biol. 2012,15;7(6):1040-8.2. Kumar P, et al. Mol Endocrinol. 2011;25(9): 1513-26.

Note: This product is for in vitro research use only