

GCP5 Polyclonal Antibody

Catalog No: #40957

Package Size: #40957-1 50ul #40957-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	GCP5 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC
Species Reactivity	Hu Ms Rt
Specificity	GCP5 Polyclonal Antibody detects endogenous levels of GCP5 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human GCP5.
Target Name	GCP5
Other Names	TUBGCP5; GCP5; KIAA1899; Gamma-tubulin complex component 5; GCP-5
Accession No.	Swiss-Prot: Q96RT8NCBI Gene ID: 114791
Uniprot	Q96RT8
GeneID	114791;
SDS-PAGE MW	118kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

Images



Western Blot analysis of 293 cells using GCP5 Polyclonal Antibody



Western Blot analysis of Jurkat A549 cells using GCP5 Polyclonal Antibody

Note: This product is for in vitro research use only