


Description

Product Name	VER-49009
Brief Description	Inhibitors
Purification	99.50%
Target Name	HSP inhibitor
Calculated MW	387.82
Formulation	C19H18ClN3O4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: CCNC(=O)C1=C(C(NN1)=C1C=C(Cl)C(O)=CC1=O)C1=CC=C(OC)C=C1
 Pathways: Cytoskeletal Signaling;
 Metabolism Receptor: HSP90β
 Boiling pt: Melting pt: Solubility: DMSO: 77 mg/mL (198.54 mM)
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Sharp SY, et al. Mol Y Ther. 2007, 6(4), 1198-1211.

Note: This product is for in vitro research use only