

## GCP6 Polyclonal Antibody

Catalog No: #40958

Package Size: #40958-1 50ul #40958-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	GCP6 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC
Species Reactivity	Hu
Specificity	GCP6 Polyclonal Antibody detects endogenous levels of GCP6 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human GCP6.
Target Name	GCP6
Other Names	TUBGCP6; GCP6; KIAA1669; Gamma-tubulin complex component 6; GCP-6
Accession No.	Swiss-Prot: Q96RT7NCBI Gene ID: 85378
Uniprot	Q96RT7
GeneID	85378;
SDS-PAGE MW	200kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

## Images



Western Blot analysis of COLO205 cells using GCP6 Polyclonal Antibody

---

Note: This product is for in vitro research use only