

## Trans-Tranilast

Catalog No: #S2270


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Trans-Tranilast
Brief Description	Inhibitors
Purification	98.00%
Target Name	RAAS antagonist
Calculated MW	327.33
Formulation	C18H17NO5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Nervous

system SMILES: COC1=CC(=C(C=C1)/C=C/C(=O)NC2=CC=CC=C2C(=O)O)OC Pathways: Endocrinology/Hormones Receptor: Angiotensin

Receptor Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Rogosnitzky, MR, et al. *AntiY Res*, 2012. 32(7): p. 2471-8. Swiderski, K, et al. *Fibrogenesis Tissue Repair*, 2014. 7(1): p. 1.

Note: This product is for in vitro research use only