


## Description

Product Name	BPR1J-097
Brief Description	Inhibitors
Purification	98.00%
Target Name	FLT inhibitor
Calculated MW	516.62
Formulation	C27H28N6O3S
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:CancerSMILES:CN1CCN(CC1)C1=CC=C(C=C1)C(=O)NC1=NNC(=C1)C1=CC(NS(=O)(=O)C2=CC=CC=C2)=CC=C1Pathways:Angiogenesis; Tyrosine Kinase/AdaptorsReceptor:FLT3Boiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:Faint Yellow PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Lin WH, et al. Br J Y. 2012 Jan 31;106(3):475-81.

Note: This product is for in vitro research use only