

Piboserod(SB-207266)

Catalog No: #S2278


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Piboserod(SB-207266)
Brief Description	Inhibitors
Purification	99.00%
Target Name	5-HT Receptor antagonist
Calculated MW	369.51
Formulation	C22H31N3O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Nervous systemSMILES:CCCCN1CCC(CNC(=O)C2=C3OCCCN3C3=CC=CC=C23)CC1Pathways:GPCR/G Protein;
NeuroscienceReceptor:5-HT4Boling pt:Melting pt:570.3CSolubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility ,
please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Kjekshus JK, et al. Eur J Heart Fail. 2009 Aug;11(8):771-8.2. Nirogi R, et al. J Pharm Pharmacol. 2013 May;65(5):704-12.

Note: This product is for in vitro research use only