

## Endoxifen

Catalog No: #S2280


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Endoxifen
Brief Description	Inhibitors
Purification	98.00%
Target Name	EGFR inhibitor; Estrogen/progestogen Receptor agonist
Calculated MW	373.49
Formulation	C25H27NO2
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Endocrine system  
 SMILES:CC/C(=C(\C1=CC=C(C=C1)O)/C2=CC=C(C=C2)OCCNC)/C3=CC=CC=C3  
 Pathways:Angiogenesis; Endocrinology/Hormones; JAK/STAT signaling; Tyrosine Kinase/AdaptorsReceptor:hERG Potassium Channel; ER; ERR  
 Boling pt:Melting pt:519.3C  
 Solubility:DMSO: 10 mM  
 Appearance:Purple Powder  
 Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.  
 Stock solution can be stored below -20 °C for several month.

## References

1. Chae YJ, et al. Eur J Pharmacol. 2015 Apr 5;752:1-7.2. Wu X, et al. Breast Y Res. 2011 Mar 10;13(2):R27.

Note: This product is for in vitro research use only