

## Motesanib (AMG-706)

Catalog No: #S2288


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Motesanib (AMG-706)
Brief Description	Inhibitors
Purification	99.00%
Target Name	VEGFR inhibitor
Calculated MW	373.46
Formulation	C22H23N5O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Cancer  
 SMILES: CC1(C)CNC2=C1C=CC(NC(=O)C1=C(NCC3=CC=NC=C3)N=CC=C1)=C2  
 Pathways: Angiogenesis; Tyrosine Kinase/Adaptors Receptor: VEGFR1; VEGFR2; VEGFR3  
 Boiling pt: Melting pt: Solubility: DMSO: 10 mM  
 Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Polverino A, et al. *Y Res*, 2006, 66(17), 8715-8721.

Note: This product is for in vitro research use only