

GDF-15 Polyclonal Antibody

Catalog No: #40960

Package Size: #40960-1 50ul #40960-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GDF-15 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	GDF-15 Polyclonal Antibody detects endogenous levels of GDF-15 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human GDF-15.
Target Name	GDF-15
Other Names	GDF15; MIC1; PDF; PLAB; PTGFB; Growth/differentiation factor 15; GDF-15; Macrophage inhibitory cytokine 1; MIC-1; NSAID-activated gene 1 protein; NAG-1; NSAID-regulated gene 1 protein; NRG-1; Placental TGF-beta; Placental bone morphogenetic
Accession No.	Swiss-Prot: Q99988NCBI Gene ID: 9518
Uniprot	Q99988
GeneID	9518;
SDS-PAGE MW	37kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

ELISA: 1/40000.

Not yet tested in other applications.

Images

Western Blot analysis of HepG2 cells using GDF-15 Polyclonal Antibody



Note: This product is for in vitro research use only