

## 8-Hydroxyquinoline-2-carboxaldehyde

Catalog No: #S2290


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	8-Hydroxyquinoline-2-carboxaldehyde
Brief Description	Inhibitors
Purification	98.48%
Target Name	Others
Calculated MW	173.17
Formulation	C10H7NO2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
SMILES: OC1=CC=CC2=C1N=C(C=O)C=C2  
Pathways: Others  
Receptor: Others  
Boiling pt: 392.2°C  
Melting pt: 97-100°C  
Solubility: DMSO: 10 mM  
Appearance: Yellow Powder  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

N/A

Note: This product is for in vitro research use only