

## 2,8-Quinolinediol

Catalog No: #S2291


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	2,8-Quinolinediol
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	161.16
Formulation	C <sub>9</sub> H <sub>7</sub> NO <sub>2</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
SMILES: OC1=CC=CC2=C1NC(=O)C=C2  
Pathways: Others  
Receptor: Others  
Boiling pt: 403.2°C  
Melting pt: 290°C  
Solubility: DMSO: 10 mM  
Appearance: Yellow Powder  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Inagami K, et al. J Biol Chem. 1965 Sep;240(9):3682-4.

Note: This product is for in vitro research use only