

## GFAP Polyclonal Antibody

Catalog No: #40961

Package Size: #40961-1 50ul #40961-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	GFAP Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu
Specificity	GFAP Polyclonal Antibody detects endogenous levels of GFAP protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human GFAP around the non-phosphorylation site of S38.
Target Name	GFAP
Other Names	GFAP; Glial fibrillary acidic protein; GFAP
Accession No.	Swiss-Prot: P14136NCBI Gene ID: 2670
Uniprot	P14136
GeneID	2670;
SDS-PAGE MW	50kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/5000.

Not yet tested in other applications.

## Images



Western Blot analysis of RAT-MUSCLE cells using GFAP Polyclonal Antibody

---

Note: This product is for in vitro research use only