


## Description

Product Name	BMS-833923
Brief Description	Inhibitors
Purification	99.47%
Target Name	Hedgehog/Smoothened antagonist
Calculated MW	473.57
Formulation	C30H27N5O
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:CancerSMILES:CNCC1=CC(NC(=O)C2=CC=C(NC3=NC4=CC=CC=C4C(=N3)C3=CC=CC=C3)C=C2)=C(C)C=C1Pathways:GPCR/G Protein; Stem CellsReceptor:SmoothenedBoiling pt:Melting pt:Solubility:DMSO: 95 mg/mL (200.6 mM)Appearance:Yellow PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Tang JY, et al. Semin Cutan Med Surg. 2011, 30(4 Suppl), S14-18.2. Zaidi AH, et al. Y Invest. 2013, 31(7), 480-489.

Note: This product is for in vitro research use only