

GHRH-R Polyclonal Antibody

Catalog No: #40962

Package Size: #40962-1 50ul #40962-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GHRH-R Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF ELISA
Species Reactivity	Hu Ms Rt
Specificity	GHRH-R Polyclonal Antibody detects endogenous levels of GHRH-R protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human GHRH-R.
Target Name	GHRH-R
Other Names	GHRHR; Growth hormone-releasing hormone receptor; GHRH receptor; Growth hormone-releasing factor receptor; GRF receptor; GRFR
Accession No.	Swiss-Prot: Q02643NCBI Gene ID: 2692
Uniprot	Q02643
GeneID	2692;
SDS-PAGE MW	47kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of L929 COLO 293T HELA
MOUSE-brain cells using GHRH-R Polyclonal Antibody

Note: This product is for in vitro research use only