

Loxiglumide

Catalog No: #S2312


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Loxiglumide
Brief Description	Inhibitors
Purification	99.00%
Target Name	cholecystokinin antagonist
Calculated MW	461.38
Formulation	C21H30Cl2N2O5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology SMILES: CCCCCN(CCCOC)C(=O)C(CCC(O)=O)NC(=O)C1=CC(Cl)=C(Cl)C=C1 Pathways: GPCR/G

Protein Receptor: CCKB Binding pt: Melting pt: Solubility: DMSO: > 10mM Appearance: Liquid Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Teysse S, et al. *Pancreas*. 1996 Nov;13(4):407-16. 2. Shiratori K, et al. *Pancreas*. 2002 Jul;25(1):e1-5.

Note: This product is for in vitro research use only