

GSK-LSD1 dihydrochloride

Catalog No: #S2315


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	GSK-LSD1 dihydrochloride
Brief Description	Inhibitors
Purification	99.50%
Target Name	Histone Demethylase inhibitor
Calculated MW	289.24
Formulation	C14H22Cl2N2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: Cl.C1C1C(NC2CCNCC2)C1C1=CC=CC=C1
 Pathways: Chromatin/Epigenetic
 Receptor: LSD1
 Boiling pt: Melting
 pt: Solubility: DMSO: 10 mM
 Appearance: Off-white Color
 Crystal Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.
 Stock solution can be stored below -20 °C for several month.

References

1. Structural Genomics Consortium.

Note: This product is for in vitro research use only