

Acotiamide hydrochloride

Catalog No: #S2319


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Acotiamide hydrochloride
Brief Description	Inhibitors
Purification	99.00%
Target Name	AChR inhibitor
Calculated MW	541.06
Formulation	C ₂₁ H ₃₇ ClN ₄ O ₈ S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism

system SMILES: O.O.O.Cl.COC1=C(OC)C=C(C(=O)NC2=NC(=CS2)C(=O)NCCN(C(C)C)C(C)C)C(O)=C1 Pathways: Neuroscience Receptor: AChE Bolin
 g pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Crystal Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Kawachi M, et al. Eur J Pharmacol. 2011 Sep; 666(1-3): 218-25.

Note: This product is for in vitro research use only