

Indacaterol

Catalog No: #S2320


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Indacaterol
Brief Description	Inhibitors
Purification	98.00%
Target Name	Adrenergic Receptor agonist
Calculated MW	392.49
Formulation	C ₂₄ H ₂₈ N ₂ O ₃
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:Respiratory systemSMILES:CCC1=C(C=C2CC(CC2=C1)NC[C@@H](C3=C4C=CC(=O)NC4=C(C=C3)O)O)CCPathways:GPCR/G
 ProteinReceptor:β-adrenoceptorBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:Yellow PowderRemark:For obtaining a higher solubility ,
 please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Battram C, et al. J Pharmacol Exp Ther, 2006. 317(2): p. 762-70.

Note: This product is for in vitro research use only