

Fimasartan(BR-A-657)

Catalog No: #S2322


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Fimasartan(BR-A-657)
Brief Description	Inhibitors
Purification	98.00%
Target Name	RAAS antagonist
Calculated MW	501.65
Formulation	C27H31N7OS
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Cardiovascular

SystemSMILES:CCCCC1=NC(C)=C(CC(=S)N(C)C)C(=O)N1CC1=CC=C(C=C1)C1=CC=CC=C1C1=NNN=N1Pathways:Endocrinology/HormonesReceptor:AGTR1Boiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:Faint Gray PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Ryu S, et al. Biol Pharm Bull. 2013;36(3):467-74.

Note: This product is for in vitro research use only