

Radotinib

Catalog No: #S2328


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Radotinib
Brief Description	Inhibitors
Purification	99.97%
Target Name	Bcr-Abl inhibitor
Calculated MW	530.5
Formulation	C27H21F3N8O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Cancer
 SMILES: CC1=CN(C=N1)C1=CC(NC(=O)C2=CC(NC3=NC=CC(=N3)C3=NC=CN=C3)=C(C)C=C2)=CC(=C1)C(F)(F)F
 Pathways: Angiogenesis; Cytoskeletal Signaling
 Receptor: Bcr-Abl
 Boiling pt: Melting pt: Solubility: DMSO: 10 mM
 Appearance: Khaki Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Heo SK, et al. PLoS One. 2015 Jun 12;10(6):e0129853.

Note: This product is for in vitro research use only