

GK1 Polyclonal Antibody

Catalog No: #40964

Package Size: #40964-1 50ul #40964-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	GK1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF
Species Reactivity	Hu Ms Rt
Specificity	GK1 Polyclonal Antibody detects endogenous levels of GK1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human GK1.
Target Name	GK1
Other Names	GK; Glycerol kinase; GK; Glycerokinase; ATP:glycerol 3-phosphotransferase
Accession No.	Swiss-Prot: P32189NCBI Gene ID: 2710
Uniprot	P32189
GeneID	2710;
SDS-PAGE MW	57kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Immunofluorescence: 1/200 - 1/1000.

Images



Western Blot analysis of Jurkat cells using GK1 Polyclonal Antibody

Note: This product is for in vitro research use only