

## Rilpivirine

Catalog No: #S2330


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Rilpivirine
Brief Description	Inhibitors
Purification	98.00%
Target Name	Reverse Transcriptase inhibitor
Calculated MW	366.42
Formulation	C22H18N6
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Immunology SMILES: CC1=CC(=CC(=C1NC2=NC(=NC=C2)NC3=CC=C(C=C3)C#N)C)/C=C/C#N Pathways: Microbiology&virology Receptor: Reverse Transcriptase Boiling pt: Melting pt: Solubility: DMSO: 73 mg/mL (199.22 mM) Appearance: Yellow Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Janssen PA, et al. J Med Chem. 2005, 48(6), 1901-1909.

Note: This product is for in vitro research use only