

Ebastine

Catalog No: #S2335


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ebastine
Brief Description	Inhibitors
Purification	98.00%
Target Name	Histamine Receptor antagonist
Calculated MW	469.66
Formulation	C32H39NO2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Nervous systemSMILES:CC(C)(C)C1=CC=C(G=C1)C(=O)CCCN1CCC(CC1)OC(C1=CC=CC=C1)C1=CC=CC=C1Pathways:GPCR/G Protein; Immunology/Inflammation; NeuroscienceReceptor:H1 receptorBoiling pt:Melting pt:Solubility:Ethanol: 11 mg/mL (23.42 mM); DMSO: 2 mg/mL (4.25 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Ko CM, et al. J Pharmacol Exp Ther, 1997, 281(1), 233-244.

Note: This product is for in vitro research use only