

Menatetrenone

Catalog No: #S2336


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Menatetrenone
Brief Description	Inhibitors
Purification	98.00%
Target Name	Vitamin inhibitor
Calculated MW	444.65
Formulation	C31H40O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Metabolism

systemSMILES:CC1=C(C(=O)C2=CC=CC=C2C1=O)C/C=C(\C)/CC/C=C(\C)/CC/C=C(\C)/CCC=C(C)CPathways:MetabolismReceptor:VK2

analogBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:Yellow PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Ishizuka M, et al. AntiY Res. 2012 Dec;32(12):5415-20.

Note: This product is for in vitro research use only